



Water-Soluble Vitamin Analysis on ZirChrom®-SAX

Dr. Bingwen Yan and Dr. Clayton McNeff
ZirChrom Separations, Inc.

Technical Bulletin #305

Traditionally the analysis of water-soluble vitamins by reversed-phase HPLC has been complicated by the lack of retention for these compounds on conventional silica C18 columns. Other analytical approaches, such as ion-pair chromatography, have also failed to yield successful and reproducible results. Here we demonstrate efficient baseline resolution of six water-soluble vitamins in six minutes using a ZirChrom®-SAX column. This method can be combined with ZirChrom's ProTain® In-Line Protein Removal System for the analysis of these compounds in biological samples.

Introduction

In this application note we focus on the HPLC analysis of Vitamin C and five B-complex vitamins; Vitamin B₁ (thiamine), Vitamin B₂ (riboflavin), Vitamin B₃ (nicotinic acid form), Vitamin B₃ (nicotinamide form), and Vitamin B₆ (pyridoxine). All of these vitamins are water-soluble.

Chromatographers oftentimes struggle in their attempts to analyze water-soluble vitamins by HPLC. Many water-soluble vitamins are very polar. Thiamin (Vitamin B₁), pyridoxine (Vitamin B₆) and ascorbic acid (Vitamin C), for example, show almost no retention on conventional C18 columns. Reversed-phase analytical methods employing ion-pair reagents have been offered as a potential solution to this problem, but these methods tend to suffer from column-to-column reproducibility problems due to the somewhat unpredictable way ion-pairing reagents interact with the silica surface and the bonded phase.

In this technical bulletin we present a unique method for the analysis of water-soluble vitamins using a ZirChrom®-SAX HPLC column. ZirChrom®-SAX, an anion exchange material, is polyethyleneimine-coated zirconia containing a substantial amount of hydrophobic cross-linker which imparts both ion-exchange and reversed-phase characteristics. The mixed mode retention characteristics of the ZirChrom®-SAX column create the unique selectivity ideal for this application (Figure 1).

For the analysis of water-soluble vitamins in serum, or other samples containing biological matrices, we recommend the addition of the ProTain® In-Line Protein Removal System; consisting of one guard holder (part# 850-00-2) and a set of three ProTain® inserts (part# PT01-0246). Please see technical bulletins #275 and #291 for further information on the use of ZirChrom's ProTain® In-Line Protein Removal System.

Experimental

Six water-soluble vitamin standards were prepared in an aqueous solution and injected on a ZirChrom®-SAX column. The separation conditions are as follows.

Column: ZirChrom®-SAX, 150 x 4.6 mm i.d.
(part number: ZR06-1546)
Mobile Phase: 50 mM Ammonium dihydrogenphosphate,
pH 4.5
Flow rate: 1.0 ml/min.
Temperature: 30 °C
Injection Vol.: 5.0 µl
Detection: UV at 254 nm

The separation is shown in Figure 1. Under these conditions the separation of Vitamin C and the four B-complex vitamins is achieved, with good peak shape and baseline resolution, in 6 minutes.




Figure 1. Analysis of Water-Soluble Vitamins.
1=Thiamine (Vit. B₁), 2=Pyridoxine (Vit. B₆),
3= Nicotinamide (form of Vit. B₃), 4=Riboflavin (Vit. B₂),
5=Nicotinic acid (form of Vit. B₃), 6=Ascorbic acid (Vit. C)

ZirChrom columns combine the high efficiency usually associated with silica columns with complete chemical and thermal stability.

ZirChrom Separations, Inc.
617 Pierce Street, Anoka, MN 55303
1-866-STABLE-1
support@zirchrom.com