



ZirChrom®


Applications of Sub-2 μ m Zirconia-PBD Columns at Elevated pH and Temperature

**Dan Nowlan, Bingwen Yan, Clayton V. McNeff and
Richard A. Henry**

**ZirChrom Separations, Inc. 617 Pierce St., Anoka, MN
55303**




Multi-Mode Behavior of Zirconia



- Zirconia substrate exhibits polar and ionic solute interaction: **especially cation-exchange.**
- With stable organic coatings, reproducible **reversed-phase** behavior can be added.
- Extreme resistance to temperature, pH and mechanical stress are unique advantages.



Addition of RP Behavior with Coated Zirconia Phases




- Retention (and selectivity) of ionic analytes modulated by pH, buffer/salt type and concentrations, and temperature.
- Retention of neutral solutes modulated by organic solvent.



C18 and Zr-PBD are Orthogonal for Basic Drugs²

C18 (RP) columns separate mainly by hydrophobic forces and Zr-PBD columns separate by a combination of ionic and hydrophobic forces



Zr-PBD and Si-C18 have very different selectivity for ionic drugs (especially in phosphate) due to the SCX ZrO₂ component.


Solutes

1. Chlordiazepoxide
2. Hydroxyzine
3. Buclizine
4. Thiothixene
5. Doxepin
6. Amitriptyline
7. Imipramine
8. Perphenazine
9. Nortriptyline
10. Desipramine
11. Thioridazine

Data provided by Sigma-Supelco



Cation Retention Observed for Zr-PBD Even in High Organic




LC Conditions: Machine-mixed 80/20 ACN/10 mM ammonium acetate pH=6.7 without pH adjustment; Flow rate, 1.0 mL/min.; Injection volume 0.1 μ L; Temperature, 35 $^{\circ}$ C; Detection at 254 nm; Columns, Zr-PBD, 50 x 4.6 mm i.d. (3 μ m particles); Silica-C18 150 x 4.6 mm i.d., (3.5 μ m particles).



Difficult Compounds for Silica Often Separate on Zirconia


Quaternary amines paraquat and diquat are retained and resolved on Zr-PS (also Zr-PBD or bare ZrO_2) due to the **cation exchange** mechanism; 50% ACN is useful to suppress or regulate retention by RP mode.

Silica-C18:
reversed-phase (silanol tailing)



column: Discovery[®] C18, 15 cm x 4.6 mm I.D., 3 μ m
mobile phase: 5% acetonitrile in 25 mM phosphate (pH 7)
flow rate: 1 mL/min.
temp.: 35 °C
det.: UV 290 nm

Zirconia-PS:
primarily ion-exchange




column: Discovery[®] Zr-PS, 7.5 cm x 4.6 mm, 3 μ m
mobile phase: 50% acetonitrile in 25 mM phosphate (pH 7)
flow rate: 3 mL/min.
temp.: 65 °C
det.: UV 290 nm

Data provided by Sigma-Supelco



Analytical Diameter Porous Zirconia Particles




1 μm 25000X

- Particles produced by a sol-gel process with 1000Å sol
- Pore diameter 250-300Å
- Density: 2.6 g/cc (2.5X silica)
- Surface area: 25 m²/g
- Particle diameters: 3 μm and **sub-2 μm**
- Totally porous (porosity: 0.45)



Idealized van Deemter Plots

$$H = A + B/v + Cv$$



Elements of the drawing provided by Sigma-Supelco



Flow Studies on 3 μ m Zr-PBD: Alkylbenzenes




Plate height based on van Deemter Equation vs linear velocity at various temperatures for retained solutes: Alkylbenzenes, Temperature: 30 °C, Mobile phase: 55/45 ACN/water, Column: ZirChrom[®]-PBD, 50 x 4.6mm, Agilent 1100/UV with micro cell (0.007" i.d. tubing).



Flow Studies on Sub-2 μm Zr-PBD: Alkylbenzenes





Plate height based on van Deemter Equation vs linear velocity for retained solutes: Alkylbenzenes, Temperature 30 °C, Mobile phase: 50/50 ACN/water (keep k in the same range as 3 μm particles), Column: 50 x 4.6mm, Agilent 1100/UV Micro Cell/0.007" i.d. tubing.



Sub-2 μm Pressure Drop at Different Temperatures*




* 3 μm particles show about half the pressure drop



ZirChrom®

Optimization and Configuration for Elevated Temperature Operation





ZirChrom®

Background Pressure Drop Across Agilent 1100 at High Flow Rate

100% H ₂ O at 30 °C		100% H ₂ O at 75 °C	
Flow (mL/min)	BP (bar)	Flow (mL/min)	BP (bar)
1	26	1	21
2	50	2	39
3	77	3	60

* Reference point: Waters Acquity (0.005" ID inlet/0.0025" ID outlet), 60/40 ACN/water, 0.5 mL/min, **background pressure = 1700 psi (113 bar)**.



Flow Studies on Sub- $2\mu\text{m}$ Zr-PBD: Factory Instrument




Plate height based on van Deemter Equation vs linear velocity for retained solutes: Alkylbenzenes, Temperature 30 °C, Mobile phase: 50/50 ACN/water (keep k' in the same range as $3\mu\text{m}$ particles), Column: 50 x 4.6mm, Agilent 1100/UV with Standard Cell and 0.007" i.d. tubing.



Flow Studies on Sub-2 μm Zr-PBD: Factory + Micro Cell Only




Plate height based on van Deemter Equation vs linear velocity for retained solutes: Alkylbenzenes, Temperature 30 °C, Mobile phase: 50/50 ACN/water (keep k' in the same range as 3 μm particles), Column: 50 x 4.6mm, Agilent 1100/UV with Micro Cell and 0.007'' i.d. tubing.



Flow Studies on Sub-2 μm Zr-PBD: Micro Cell + Optimized Tubing




Plate height based on van Deemter Equation vs linear velocity for retained solutes: Alkylbenzenes, Temperature 30 °C, Mobile phase: 50/50 ACN/water (keep k' in the same range as 3 μm particles), Column: 50 x 4.6mm, Agilent 1100/UV with Micro Cell and optimized 0.005'' i.d. tubing.



Flow Studies on Sub-2 μm Zr-PBD: Heat Exchanger + Fitting + μCell




Plate height based on van Deemter Equation vs linear velocity for retained solutes: Alkylbenzenes, Temperature 30 °C, Mobile phase: 50/50 ACN/water (keep k' in the same range as 3 μm particles), Column: 50 x 4.6mm, Agilent 1100/UV with Micro Cell, high pressure fitting and passing through heat exchanger.



Flow Studies on Sub-2 μm Zr-PBD: Factory Instrument at Ambient




Plate height vs linear velocity for retained solutes: Alkylbenzenes, Temperature 30 °C, Mobile phase: 50/50 ACN/water (keep k' in the same range as 3 μm particles), Column: 50 x 4.6mm, Agilent 1100/UV with Standard Cell and 0.007" i.d. tubing.



Flow Studies on Sub-2 μm Zr-PBD: Micro Cell




Plate height vs linear velocity for retained solutes: Alkylbenzenes, Temperature 30 °C, Mobile phase: 50/50 ACN/water (keep k' in the same range as 3 μm particles), Column: 50 x 4.6mm, **Agilent 1100/UV with Micro Cell and 0.007" i.d. tubing.**



Flow Studies on Sub-2 μ m Zr-PBD: Micro Cell + Tubing




Plate height vs linear velocity for retained solutes: Alkylbenzenes, Temperature 30 °C, Mobile phase: 50/50 ACN/water (keep k' in the same range as 3 μ m particles), Column: 50 x 4.6mm, **Agilent 1100/UV with Micro Cell and 0.005" i.d. tubing.**



Flow Studies on Sub-2 μ m Zr-PBD: Heat Exchanger + Fitting + μ Cell




Plate height vs linear velocity for retained solutes: Alkylbenzenes, Temperature 30 °C, Mobile phase: 50/50 ACN/water (keep k' in the same range as 3 μ m particles), Column: 50 x 4.6mm, **Agilent 1100/UV with Micro Cell, high pressure fitting and passing through heat exchanger.**



Flow Studies on 3 μ m Zr-PBD: Factory Instrument




Plate height based on van Deemter Equation vs linear velocity at various temperatures for retained solutes: Alkylbenzenes, Temperature: 30 °C, Mobile phase: **50/50 ACN/water**, Column: ZirChrom[®]-PBD, 50 x 4.6mm, Agilent 1100/UV standard cell (0.007" i.d. tubing).



Flow Studies on 3 μ m Zr-PBD: Factory + Micro Cell





Plate height based on van Deemter Equation vs linear velocity at various temperatures for retained solutes: Alkylbenzenes, Temperature: 30 °C, Mobile phase: 50/50 ACN/water, Column: ZirChrom®-PBD, 50 x 4.6mm, Agilent 1100/UV micro cell (0.007" i.d. tubing).



ZirChrom®

Alkylbenzylamine Separation on **sub-2 μ m Zr-PBD, 50 °C**




LC Conditions: Column: ZirChrom®-PBD, 50 x 4.6 mm i.d., sub-2 μ m Mobile Phase: 21/79 ACN/20 mM K₃PO₄ at pH=12; Flow rate: 1.5 mL/min; Temperature: 50 °C; Injection Vol.: 3.0 μ L; Detection: UV at 254 nm



ZirChrom®

Antihistamine Separation on **sub-2 μ m Zr-PBD, 80 °C**



LC Conditions: Column: ZirChrom®-PBD, 50 x 4.6 mm i.d., sub-2 μ m Mobile Phase: 28/72 ACN/50 mM TMA-OH at pH=12.2; Flow rate: 2.5 mL/min; Temperature: 80 °C; Injection Vol.: 2.0 μ L; Detection: UV at 254 nm



ZirChrom®

β -blockers on ZirChrom®-PBD sub-2 μ m, High Temp

Analytes

- 1=Labetalol
- 2=Atenolol
- 3=Acebutolol
- 4=Metoprolol
- 5=Oxprenolol
- 6=Lidocaine
- 7=Quinidine
- 8=Alprenolol
- 9=Propranolol

ZirChrom®-PBD


50mm x 4.6mm, sub-2 μ m

22/78 ACN/20mM K₃PO₄ at pH=12

F=2.5 mL/min

UV=254nm

T=75 °C





ZirChrom®

β -blockers on ZirChrom®-PBD sub-2 μ m, High Temp, Faster sampling

Analytes

- 1=Labetalol
- 2=Atenolol
- 3=Acebutolol
- 4=Metoprolol
- 5=Oxprenolol
- 6=Lidocaine
- 7=Quinidine
- 8=Alprenolol
- 9=Propranolol

ZirChrom®-PBD


50mm x 4.6mm, sub-2 μ m

21/79 ACN/20mM K₃PO₄ at pH=12

F=2.5 mL/min

UV=254nm

T=75 °C, faster sampling rate





Plans for Further Development

- Extend the range of ultra-high speed applications using sub-2 μ m Zr-PBD, especially at high pH and temperature (“extreme conditions for silica”); develop generic conditions for LC-MS.
- Develop sub-2 μ m Zr-CARB and compare performance to Zr-PBD under ambient and extreme temperature conditions.
- Study additional advantages of optimizing the IBW of an Agilent Model 1100 HPLC instrument using a high performance (Model 1200) heat exchanger.



ZirChrom®

Conclusions

- Multi-mode HPLC columns have become popular for difficult applications where compounds have ionic character and vary widely in chemical nature. Several ZirChrom® phases are ideal and popular for multi-mode applications and are stable over a much wider range of pH and temperature than any silica-based phase.
- Zirconia 3 μ m HPLC columns are currently available in a wide range of stable coatings and produce efficiencies in excess of 100,000 N/M.
- New sub-2 μ m zirconia UHPLC columns with very high efficiency in excess of 200,000 N/M with a PBD polymer-coated phase permit higher speed separations with shorter residence time at elevated temperatures.



References and Acknowledgements

1. **J. A. Blackwell and P. W. Carr, "Development of an Eluotropic Series for the Chromatography of Lewis Bases on Zirconium Oxide," Anal. Chem. 64, 863-73 (1992).**
2. **R. A. Henry, H. K. Brandes, D. S. Bell and C. T. Santasania, 30th Annual HPLC Meeting, Oral Presentation, 2006, San Francisco, CA.**
3. **R. A. Henry and H. K. Brandes, Eastern Analytical Symposium, Oral Presentation, 2006, Somerset, NJ.**
4. **B. Yan. C. V. McNeff. R. A. Henry and D. Nowlan, Eastern Analytical Symposium, Poster, 2008, Somerset, NJ.**
5. **B. Yan. C. V. McNeff. R. A. Henry and D. Nowlan, Pittsburgh Analytical Conference, 2009, Oral Paper, Chicago, IL..**
6. **B. Yan. C. V. McNeff. R. A. Henry and D. Nowlan, Horvath Award Symposium, 2009, Oral Paper, Hartford, CT.**

The assistance of Supelco Division of Sigma-Aldrich is gratefully appreciated, including the use of a high-pressure column fitting.