



Application of Zirconia Phases for Structurally Similar Compounds and other Difficult HPLC Separations

Kelly S. Johnson¹, Merlin Bicking², and Richard A. Henry³

¹ZirChrom Separations, 617 Pierce St, Anoka, MN 55303


²ACCTA, Inc., 3534 Jessie Ct, Saint Paul, MN 55125

³Consultant, Supelco, Div. of Sigma Aldrich

Why Zirconia?


- Bigger Method Development Toolbox
- Buffer Type and Concentration can "Tune" Selectivity
- Stable Over a Wide Range of pH (1-14)
- Stable at Temperatures Up to 200 °C
- Selectivity
- Unique Multimodal Selectivity

Zirconia Surface Chemistry



Ligand Exchange Interaction

High pH Stability Comparison*




Exposure Conditions: Mobile phase, ACN/50 mM Potassium phosphate buffer at indicated pH; Temperature, 30 °C.
LC Conditions: Mobile phase, ACN (or THF)/50 mM Potassium phosphate buffer at indicated pH; Flow Rate, 1.0 mL/min.; Temperature, 30 °C; Injection Volume, 5 µL; Detection, 254 nm.

* Column names are the trademarks of their respective manufacturers.

Elevated Temperature – Enables Fast Separation Without Switching Columns or Loss of Resolution

LC Conditions: Mobile Phase, 29/71 ACN/50mM Tetramethylammonium hydroxide, pH 12.2; Flow Rate, 1.35 mL/min.; Injection volume, 0.5 µL; 254 nm detection; Column Temperature, 21 °C; Pressure drop = 195 bar; Solutes: 1=Doxylamine, 2=Methapyrilene, 3=Chlorpheniramine, 4=Triprolidine, 5=Meclizine 100 x 4.6 ZirChrom-PBD



LC Conditions: Mobile Phase, 26.5/73.5 ACN/50mM Tetramethylammonium hydroxide, pH 12.2; Flow Rate, 3.00 mL/min.; Injection volume, 0.2 µL; 254 nm detection; Column Temperature, 80 °C; Pressure drop = 195 bar; Solutes: 1=Doxylamine, 2=Methapyrilene, 3=Chlorpheniramine, 4=Triprolidine, 5=Meclizine 100 x 4.6 ZirChrom-PBD




Choosing Buffer Type


Interaction Strength	Lewis Base Additive (A)	Retention
Strongest	Hydroxide	More Retention
	Phosphate	A stronger interaction strength means a greater charge on the surface, producing more retention for cations.
	Fluoride	Anions would be affected in the opposite way, and neutrals are not affected.
	Citrate	Less Retention
	Sulfate	
	Acetate	
	Formate	
Weakest	Chloride	

Small Lewis bases with higher electron density and lower polarizability interact more strongly with zirconia.


Effect of Buffer Type & Concentration on Selectivity



Effect of pH on Selectivity




ZirChrom®-SAX



ZirChrom®-SAX is a zirconia particle coated with polyethyleneimine polymer and then crosslinked using a hydrophobic crosslinker. The crosslinker modifies the phase with a bit of reversed phase character, giving the ZirChrom-SAX phase a unique, multimodal selectivity.

Water Soluble Vitamins on ZirChrom®-SAX




1=Thiamine (Vit. B₁), 2=Pyridoxine (Vit. B₆), 3=Nicotinamide (form of Vit. B₃), 4=Riboflavin (Vit. B₂), 5=Nicotinic acid (form of Vit. B₃), 6=Ascorbic acid (Vit. C)

Six water soluble vitamins are easily separated under isocratic conditions.

Column: ZirChrom®-SAX, 150 x 4.6 mm i.d., 3 µm
 Mobile Phase: 50 mM Ammoniumdihydrogenphosphate, pH 4.5
 Flow rate: 1.0 ml/min.
 Temperature: 30 °C
 Injection Vol.: 5.0 µl
 Detection: UV at 254 nm

“Green” Analysis of Energy Drink




Columns: Leading C18 column, ZirChrom®-SAX, 100 x 4.6 mm i.d., 3 µm
 Mobile Phase: 10 mM NH₄H₂PO₄ + 5 mM NH₄CO₃
 Flow rate: 1.5 ml/min
 Temperature: 50 °C
 Injection Vol.: 5 µL
 Detection: UV at 230 nm


Traditional C18 phases do not readily separate benzoate and sorbate (preservatives in energy drinks). The unique multimodal selectivity of ZirChrom-SAX can take advantage of differences between the ionized forms, producing excellent selectivity and retention, resulting in 6 minute baseline resolution of caffeine (C), benzoate (B) and sorbate (S).

“Green” Mobile Phase

Chemical Structure of ZirChrom®-CARB




Comparison of Selectivity of ODS, ZirChrom®-PBD &-CARB



- N-benzyl formamide
- Benzylalcohol
- Phenol
- 3-phenyl propanol
- p-chlorophenol
- Acetophenone
- Benzonitrile
- Nitrobenzene (NB)
- methyl benzoate
- Anisole
- Benzene
- p-chlorotoluene
- p-nitrobenzyl chloride
- Toluene
- Benzophenone
- Bromobenzene
- Naphthalene
- Ethylbenzene
- p-xylene
- p-dichlorobenzene
- Propylbenzene
- n-butylbenzene

Structurally Similar Sulfated Steroids on ZirChrom®-CARB



Column: 4.6 mm x 100 mm, ZirChrom-CARB
 Mobile Phase: Gradient elution from 55/5/40 to 90/5/5 A/B/C from 0 to 4.5 minutes.
 A: Acetonitrile
 B: Tetrahydrofuran
 C: 25 mM Ammonium fluoride, 10 mM Ammonium acetate, pH 5.6
 Temperature: 80 °C
 Injection Vol.: 10 µl
 Flow rate: 3.0 ml/min.
 Pressure Drop: 195 bar
 Detection: UV at 270 nm

Conclusions

- Zirconia-based stationary phases have very different selectivity from silica-based phases.
- ZirChrom®-SAX has a unique multimodal selectivity – ionic and hydrophobic.
- Zirconia phases have selectivity that can be tuned by the addition of different buffers.
- Zirconia reversed phases, such as ZirChrom®-CARB offer excellent chemical and thermal stability:
- Faster separations at high temperature.
- Better selectivity for structurally similar compounds.

Acknowledgements

The authors extend their appreciation to the following sponsors who supplied resources for this research.

